



**ANTICON**<sup>®</sup>  
BY CONTEC

When it has to be clean, it has to be Anticon<sup>®</sup> by Contec



**Anticon<sup>®</sup>  
Sterile Gold<sup>™</sup>**

*Assured Sterility via Gamma Irradiation*

*Sealed Edge for Low Particles and Fibers*

*Superior Value for General Cleaning*



# Anticon<sup>®</sup> Sterile Gold<sup>™</sup>

**100% sterile polyester double-knit interlock wipe with sealed edges**

**DESCRIPTION**

Contec's Anticon Sterile Gold is made from Anticon Gold HeavyWeight wipes. They are made from 100% continuous filament polyester with a double-knit no-run interlock construction. Anticon Sterile Gold uses patented textile technology to achieve superior sorption characteristics for critical applications. The proprietary sealed edges of these wipers allow the use of 100% of the wipe for cleaning without the worry of harming delicate surfaces. Anticon Sterile Gold wipes are ideal for use in the most critical aseptic environments. Anticon Sterile Gold is conveniently packaged in easy tear bags containing 20 wipes per package. Anticon Sterile Gold is sterilized via gamma irradiation, and a certificate of irradiation is provided with each shipment.

**FEATURES & BENEFITS**

- Gamma irradiated to a Sterility Assurance Level of 10<sup>-6</sup> SAL which eliminates the introduction of biological contamination
- Sealed edges to reduce particle and fiber generation
- High sorptive capacity and fast adsorption rate for superior spill control
- Manufactured from 100% continuous filament polyester and specially processed for superior cleanliness and sorbency
- Cleanroom laundered
- Tested in accordance with the IEST Recommend Practices for the Evaluation of Wiping Materials for Critical Environments
- Ideal for applying and removing disinfectants and cleaning solutions
- Packaged in low count easy tear bags to reduce waste
- Designed for use in aseptic environments

**ORDERING INFORMATION**

Part #	Description	Packaging
492226-942	12" x 12" (31 x 31cm) Sterile Gold wipes	20 wipes/IB, 5 IB/OB, 4 OB/case

KEY: Bulk= Advantage Pack    IB= Inner Bag    OB= Outer Bag    Blue= Sterile Part Numbers

Anticon<sup>®</sup> Cleanroom Products by Contec, Inc.  
525 Locust Grove  
Spartanburg, SC 29303

Customer Service  
1-800-762-3472  
+1-864-503-8333  
anticon@contecinc.com

 **CONTEC**<sup>®</sup>  
www.anticonwipers.com  
www.contecinc.com



**ANTICON<sup>®</sup>**  
BY CONTEC

**Anticon<sup>®</sup>**  
**Sterile Gold<sup>™</sup>**

<b>TYPICAL TEST DATA</b>			
<b>Basis Weight</b>		(gm/m <sup>2</sup> )	140
<b>Liquid Particles</b>			
Biaxial Shake Test (IEST-RP-CC-004.2 Sec 5.2)	≥ 0.5 μm	(million particles/m <sup>2</sup> )	7.74
	≥ 5.0 μm	(million particles/m <sup>2</sup> )	0.11
<b>Fiber Analysis</b>			
ASTM F311-F312	> 100 μm	(fibers/m <sup>2</sup> )	47
<b>NonVolatile Residue</b>			
Short Term Extraction (IEST-RP-CC-004.2 Sec 6.1.2)	DI Water	(grams/m <sup>2</sup> )	0.004
	IPA	(grams/m <sup>2</sup> )	0.030
	Acetone	(grams/m <sup>2</sup> )	0.066
<b>Sorbency</b>			
Rate & Capacity (IEST-RP-CC-004.2 Sec 7.1 & 7.2)	Rate	(seconds)	< 1
	Capacity	(ml/m <sup>2</sup> )	573
<b>Specific Ion Content</b>			
Short Term Extraction (IEST-RP-CC-004.2 Sec 6.1.2)	Na	(μg/kg (ppb))	18
	Li	(μg/kg (ppb))	< 1
	NH4	(μg/kg (ppb))	19
	K	(μg/kg (ppb))	84
	Mg	(μg/kg (ppb))	< 1
	Ca	(μg/kg (ppb))	2
	Fl	(μg/kg (ppb))	27
	Cl	(μg/kg (ppb))	10

**APPLICABILITY:** This information is applicable to the following part numbers as well as to any product containing the same materials and produced under the same conditions, regardless of product size or packaging configuration. Please contact a Contec sales representative for more details.  
492226-942

**NOTES:** This test data represents typical values and should not be used as product specifications. Valid product comparison may only be obtained through side-by-side testing in the same test facility, under similar conditions.