



When it has to be clean, it has to be Anticon® by Contec



Anticon® Millimoist® Gold

Uniform Saturation Level

Sealed Edge for Low Particles and Fibers

Eliminates Solvent Mixing and Storage

*Ask about Anticon® Millimoist Gold with
Particle Attraction Technology (P.A.T.)*

Anticon® Millimoist® Gold

Cleanliness and convenience in a premoistened wipe

DESCRIPTION

Anticon Millimoist Gold is a manufactured using our ultra-clean Anticon Gold wipes in conjunction with semiconductor grade solvents. All wipes are moistened with the optimum amount of cleaning solution to deliver absolute cleaning power, not excess solution. Millimoist Gold polyester premoistened wipes are delivered with a patented low particle generation technology. Use Millimoist with the knowledge and security that there are no undo particle burdens being added to the critical environment.

FEATURES

- 100% continuous filament polyester, chemically modified for superior cleanliness and sorbency
- Tested in accordance with the IEST Recommend Practices for the Evaluation of Wiping Materials for Critical Environments
- Consistent lot to lot performance
- Produced using the proper grade of solvent for the user's environment

BENEFITS

- Less contamination due to lowest sub-micron particles and fiber counts
- Eliminates variation in cleaning
- Eliminates storing and mixing solvents
- Reduces solvent consumption

APPLICATIONS

- Designed for use in ISO Class 1-4 environments
- Cleanroom wall wipe downs
- Superior chamber cleaning wipe
- Pre-gown wipe down of pagers, phones, pens, etc.

ORDERING INFORMATION

Part #	Description	Packaging
492239-101	9" x 9" (23 x 23cm) Gold Sorb wipes, 5% IPA/95% DI Water, bulk	50 wipes/pouch; 16 pouches/case
495357-660	9" x 9" (23 x 23cm) Gold StandardWeight wipes, 6% IPA/94% DI Water, bulk	50 wipes/pouch; 16 pouches/case
495357-070	9" x 9" (23 x 23cm) Gold StandardWeight wipes, 70% IPA/30% DI Water, bulk	50 wipes/pouch; 16 pouches/case
495357-680	9" x 9" (23 x 23cm) Gold StandardWeight wipes, 100% IPA, bulk	50 wipes/pouch; 16 pouches/case

KEY: Bulk= Advantage Pack



ANTICON[®]
BY CONTEC

Anticon[®]
Millimoist[®] Gold

TYPICAL TEST DATA			495357	492239
Basis Weight		(gm/m ²)	120	140
Liquid Particles				
Biaxial Shake Test (IEST-RP-CC-004.2 Sec 5.2)	≥ 0.5 μm	(million particles/m ²)	7.44	6.94
	≥ 5.0 μm	(million particles/m ²)	0.13	0.15
Fiber Analysis				
ASTM F311-F312	> 100 μm	(fibers/m ²)	48	79
NonVolatile Residue				
Short Term Extraction (IEST-RP-CC-004.2 Sec 6.1.2)	DI Water	(grams/m ²)	0.009	0.005
	IPA	(grams/m ²)	0.036	0.021
Specific Ion Content				
Short Term Extraction (IEST-RP-CC-004.2 Sec 6.1.2)	Na	(μg/kg (ppb))	28	22
	K	(μg/kg (ppb))	72	77
	Mg	(μg/kg (ppb))	<1	<1
	Ca	(μg/kg (ppb))	28	<1
	Cl	(μg/kg (ppb))	9	12

APPLICABILITY: This information is applicable to the following part numbers as well as to any product containing the same materials and produced under the same conditions, regardless of product size or packaging configuration. Please contact a Contec sales representative for more details.
492239-101, 495357-660, 495357-680, 495359-070

NOTES: This test data represents typical values and should not be used as product specifications. Valid product comparison may only be obtained through side-by-side testing in the same test facility, under similar conditions.